

# SONICA® 5200 S3

Product name	Code
SONICA 5200 M S3 *	090.003.0030
SONICA 5200 MH S3 *	090.003.0031
SONICA 5200 ETH S3 *	090.003.0032
SONICA 5200 EP S3 *	090.003.0033
Ultrasonic Cleaner 21 Litres	
Accessories	
Stainless steel wire basket Internal dim. 460x260x90mm	090.004.0009

	Internal dim. 460x260x90mm	
SONICA 5200	Stainless steel lid 5200/5300	090.004.0025
	SONICA 5200	090.004.0047

for two glass or plastic beaker
Plastic beaker 250 ml 090.004.0032

Further accessories on request

Stainless steel beaker holder



\* Stainless steel tank made in AISI 304 thickness 08/10 with round tank corners for a better removing of dirt residuals after cleaning process. The tank is conical and realized with a press machine

# **General Characteristics**

Glass beaker 400 ml

Ultrasonic Cleaning is a safe and efficient modern procedure, which ensures perfect cleaning within the shortest time possible. It has many advantages over ordinary cleaning methods. Ultrasonics can remove even the most tenacious deposits from parts which would be difficult to reach by hand, it also reduces the risk of injury or infection which may occur during manual cleaning.

# Notes for use as "Medical Devices"

**Denomination**: Ultrasonic Cleaning Equipment EMDN: Z12011302

Field of application: Ultrasonic Cleaning Equipment for surgical and dental instruments

090.004.0036

090.004.0033

Classification: Medical Devices Regulation Class I rule 13, Active MD, Non Invasive MD, Non Implantable

MD

Product MDR compliant Regulation (EU) 2017/745

### SOLTEC S.r.l.

Uffici Commerciali e Sede Legale: Via G. Röntgen 16 - 20136 Milano
Tel. +39 02 58 308 378 - Fax +39 02 58 308 595 - www.soltec.it - info@soltec.it
Produzione e Magazzino: Via G. Castelbarco, 17 - 20136 Milano - Tel. +39 02 58 324 131

Cod. fisc. e P.IVA: IT11127210158

N. Mecc.: MI 221750 - C.C.I.A.A. Milano N. 1441548

lscr. Trib. Mi Reg. Soc. N. 343391/8444/41

Cap. Soc. Euro 110.500 i.v. - Iscr. Registro AEE Nr. IT08020000000435





# **Ultrasonic Cleaning Systems** for Medical and Industry

Rev.4-04-2024

#### **Technical Data**

Technical specification subject to change

Main line voltage	230/240V~50/60Hz
Main input power/consumption	400 W
Main input power with heating	1400 W
Peak HF Ultrasound Power	800 W
US frequency	39 kHz ± 1 kHz with SWEEP SYSTEM Technology
Liquid drain	1/2"
Transducers	8 (2 Piezoelectric ceramics each)
Tank capacity [L]	21
Weight [KG]	14.5
Installation category	Class II according to EN 61010-1
Environmental conditions	Temperature Rel. Hum. 80% up to 31°C with linear
	decrease up to 50% at 40°C
ID annuals	ID 00

IP grade **CE-compliant** 

**IP 22** EMC (EN 61326-1)

LV (EN 61010-1 EN 61010-2-040)

2017/745/UE Medical Device Regulation Risk Class I according to the rule 13 to the MDR

Overheat switch: Automatic bimetallic Thermostat operating at 140°C fixed on the heating element


## Approximately Internal Tank Dimensions

# **External Appliance Dimensions**

W/D/H (mm)	500 / 300 / 150	W/D/H (mm)	600 / 330 / 425
Max filling volume tank	21 L	Max operating tank volume	18 L

#### **Constructive Materials**

Housing Stainless steel AISI 304 Tank Stainless steel AISI 304



# Control panel available series

M SONICA

Sweep System generator

MH

SONICA

**ETH** 

**EP** 

M serie: The simplest version in the MH serie: This version is equipped SONICA range is equipped with a timer which can be set for 15 minutes. It can also be set to manual operation up to 1 hour of continuous. No heating system

with a timer which can be set for 15 minutes. It can also be set to manual operation up to 1 hour of continuous. The thermostat for the heater is set

Sweep System generator



ETH serie: Cleaning times and EP serie: The most technologically temperature controlled electronically by a microprocessor.

minutes. Double washing times 10-20-30

minutes by pressing twice the minutes key

Sweep System generator



advanced of the SONICA range: Sweep System generator; modern Ultrasonic cleaning times 5-10-15 digital display control panel window; degas time setting; cleaning times from 0 to 99 minutes; tank liquid temperature setting up to 70°C; memory for 9 cleaning programmes; jet program function to start quickly and easily cleaning cycle.

SOLTEC Srl - Via G.Rontgen 16 - 20136 Milano - Italy- Tel +39 0258308378 - Fax +39 0258308595 info@soltec.it https://www.soltec.it